

**Quantachrome NovaWin - Data Acquisition and Reduction
for NOVA instruments
©1994-2010, Quantachrome Instruments
version 11.0**



Analysis

Operator: UEL - QUIMICA
Sample ID: LECA 28G
Sample Desc: LECA 28G
Sample weight: 0.1307 g
Outgas Time: 3.0 hrs
Analysis gas: Nitrogen
Press. Tolerance: 0.100/0.100 (ads/des)
Analysis Time: 155.9 min
Cell ID: 0

Date: 2022/11/10

Filename: C:\Users\236795\Documents\Odair\LECA 28G\LECA 28G.qps
Comment:
Sample Volume: 0 cc
Outgas Temp: 300.0 C
Bath Temp: 77.3 K
Equil time: 60/60 sec (ads/des)
End of run: 2022/11/10 14:26:41

Report



Operator: UEL - QUIMICA Date: 11/18/2022
C:\Users\236795\Documents\Odair\LECA 28G\LECA 28G.qps

Equil timeout: 240/240 sec (ads/des)
Instrument: Nova Station A
F/W version: 0.00

Multi-Point BET Plot

Data Reduction Parameters

Adsorbate	Nitrogen	Temperature	77.350K	Liquid Density:	0.808 g/cc
	Molec. Wt.: 28.013 g	Cross Section:	16.200 Å²		



BET summary

Slope = 120.160
Intercept = 1.416e+00
Correlation coefficient, r = 0.999909
C constant = 85.874
Surface Area = 28.645 m²/g