

Fisons- EA1108
Company name: CE Instruments
Software: EAGER 200
Operator ID: LV
Method in use: CHNS2401

Data de análise: 21/02/2024

Código_CA	Massa (mg)	%N	%C	%H	%S
A12240201	1,273	2,61	47,09	6,59	0
A12240202	1,220	4,28	43,33	6,69	0
A12240203	1,053	0,92	35,69	6,38	0
A12240204	1,059	0	38,54	6,68	0
A12240205	1,720	4,34	39,65	3,01	0
A12240206	2,126	2,025	33,54	3,19	0
A12240207	1,357	2,68	47,16	3,51	0
A12240208	2,009	3,19	47,23	3,01	0
A12240209	1,507	3,68	44,32	2,68	0
A12240210	1,370	2,50	46,73	3,58	0
A12240211	0,872	3,57	56,87	6,19	0
A12240212	1,628	0	51,71	7,56	0

Padrão/Controle	Massa (mg)	%N Cert. Exp.		%C Cert. Exp.		%H Cert. Exp.		%S Cert. Exp.	
Metionina	1,297	9,42	9,45	40,26	40,31	7,48	7,52	21,50	21,49
BBOT	1,656	6,48	6,51	72,51	72,54	6,08	6,11	7,42	7,46
Cistina	1,895	11,64	11,67	30,00	30,00	5,02	5,02	26,67	26,68
Sulfanilamida	1,753	16,29	16,31	41,77	41,79	4,65	4,64	18,67	18,70
Acetanilida	1,675	10,36	10,38	71,02	71,13	6,71	6,69	0	0

Cert.: especificação no certificado

Exp.: resultado experimental

Micro balança: Sartorius Micro (XM-1000 P)
Calibrada em: 24/01/2024 (Evagon - Jundiaí/SP)

Amostra	Código_CA
Leca1	A12240201
Leca2	A12240202
Leca3	A12240203
Leca4	A12240204
Leca6	A12240205
Leca7	A12240206
Leca13G	A12240207
Leca15G	A12240208
Leca28G	A12240209
Leca29G	A12240210
CHSTH-H3PO4	A12240211
bagço	A12240212